5.2 Electrometric indication methods: In the meantime indication methods using two polarized (polarizable) electrodes have established themselves. Double Pt-wire or double Pt-ring electrodes are commonly used.

5.2.1 Biamperometric indication (Upol or dead stop)

In this case a constant voltage of max. 500 mV is applied to the electrodes and the resulting current is measured. This method was first used mainly for aqueous iodometric titrations. The following reactions occur at the electrodes:

The relationship between the voltage, current and concentration of the depolarizer (iodine or triiodide) is similar to that found in polarography, i.e. initially only low voltages are required to produce a flow of current. Then the voltage increases steeply – as a function of the iodine concentration – and, when the limiting diffusion current has been achieved, the current increases only slightly. If the voltage is then further increased decomposition of either the solvent or conductive salt could occur. This means that a current flows for as long as iodine is present. When the iodine has been consumed (e.g. by titration with thiosulfate) the current falls abruptly to 0 (which is where the name «dead stop» comes from), the resistance between the electrodes increases strongly and the electrodes are again polarized. In KF titration the reverse occurs. If excess H₂O is present only a minimal current flows; this increases to a few μ A when excess iodine is present. The «water

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loss» resulting from electrochemical conversion can be ignored. At 50 μA for 30 s it is theoretically approx. 1 μg H_2O . The illustration shows an example of a biamperometrically indicated titration curve:

5.2.2 Bivoltametric indication (Ipol)

In this case a small direct or alternating current is applied between the electrodes (AC is much more sensitive than DC, which is why AC is usually only used for coulometric titrations) and the resulting voltage is measured. The shape of the titration curve is similar to that obtained by biamperometry; however, the jumps are usually more marked and larger. If the sample solution and titrant are electrochemically active then peak-shaped or Vshaped titration curves are obtained. If only one partner is electrochemically active (the KFR with excess iodine) then L-shaped curves are obtained. Examples:

The resulting voltages for excess iodine or water depend, of course, on the applied current (Ipol). Three examples (double Pt-wire electrode): Ipol = 50 μA , U approx. 100 mV for excess iodine, approx. 680 mV for excess H_2O Ipol = 10 μA , U approx. 5 mV for excess iodine, approx. 580 mV for excess H_2O Ipol = 1 μA , U approx. 0 mV for excess iodine, approx. 350 mV for excess H_2O

The titration endpoint is determined by entering a switch-off voltage; reagent addition ceases when this limit is undercut. This «ideal endpoint» is best found experimentally. If too low a

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voltage is selected then the amount of excess iodine required to achieve the switch-off point will be too large. If it is too high then the titration will not start because no free iodine is required for this voltage to be achieved (decomposition of the solvent or electrolyte). In all these indication methods there is an overtitation (presence of excess iodine). In comparison to the visual methods, the electrometric methods have the advantage that the titration is always carried out to the same (slight) excess of iodine and in this way a better reproducibility and accuracy can be achieved.

6.2 Coulometry: In contrast to volumetric KF titration, in KF coulometry no titrant is added but the required iodine is generated directly in the iodide-containing electrolyte by electrochemical means. There is a strictly quantitative relationship between the amount of electric charge and the amount of iodine generated. This relationship allows the extremely precise addition of the iodine. Metrohm KF Coulometers work with variable current strengths and pulse lengths for generating iodine at the anode. Higher current strengths have not proven themselves in practice – side reactions can occur and heat can be produced. Working with variable pulse lengths allows precise «iodine addition» even in the region of the titration endpoint. The instrument measures the time and current flow that is required to reach the titration endpoint. The product of time x current is directly proportional to the

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amount of iodine generated and therefore the amount of water determined. Coulometric water determination is primarily used for the determination of small amounts of water. Larger amounts of water require a lot of time and/or may exceed the water capacity of the KF reagent, which could lead to incorrect results. Metrohm KF Coulometers work in a determination range of 10 µg...200 mg H₂O with a resolution of 0.1 µg H₂O. The iodine is produced by a generator electrode. There are two different types of generator electrode:

A) Diaphragm-less generator electrode

The diaphragm-less generator electrode is easy to handle and easy to clean. Only a single reagent is required. However, care must be taken that only those reagents that are intended for diaphragm-less generator electrodes are used. As no moisture can adhere to the diaphragm there is a further advantage – the electrode is rapidly ready for use. There is no cathode compartment in the diaphragm-less generator electrode, as reduction takes place on a platinum surface.

B) Generator electrode with diaphragm A generator electrode with diaphragm should be used whenever your samples contain ketones and aldehydes, as the special reagents used for them are only available for generator electrodes with diaphragms. An electrode with diaphragm should also be used for reagents with a low conductivity (the addition of chloroform, e.g., reduces the conductivity) and when measurements

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are to be made in the lower **trace (vết, chút ít, số lượng thấp)** range. Reagents for coulometric water determination consist of an anolyte that is filled into the titration vessel and a catholyte that is filled into the cathode compartment. Both liquids should have about the same level in order to prevent pressure compensation from forcing the anolyte into the cathode compartment and vice versa. The endpoint recognition is by a double platinum electrode and is largely similar to that in volumetric Karl Fischer titration. As the coulometric KF method is an absolute method this means that no titer has to be determined. It must only be ensured that the reaction that generates the iodine takes place to 100%. Although this is the case with the reagents that are currently available, we still recommend that the whole system is checked from time to time with the aid of certified water standards with a water content of 1.0 ± 0.03 mg/g or 0.1 ± 0.005 mg/g. Handling the standards is described in 6.1. In coulometry, cells made completely of glass are used as titration cells, as cells with a plastic cover have a certain permeability to water vapor. This does not have much influence on volumetric KF titrations, as these are normally used to determine higher water contents than in KF coulometry. Care must also be taken that the cell is as impervious as possible.

6.3 The Karl Fischer oven Many substances only release their water slowly or at elevated temperatures. This means that they are not suitable

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for direct Karl Fischer titration. A further problem is the low solubility of certain samples in alcohols. In such cases traditional methods recommend complicated sample preparation procedures or the use of solubility promoters, which represent a health hazard. Other substances react with the KF reagents to release water or consume iodine; this falsifies the results. The above-mentioned problems can be avoided by using a heating method.

The substance under investigation is heated in an oven located upstream from the cell and the released water is transferred by a flow of dry carrier gas to the titration cell where it is determined by a Karl Fischer titration. As only the water enters the titration cell and the sample itself does not come into contact with the KF reagent this means that unwanted side reactions and matrix effects are ruled out. Work is carried out with either a conventional tube oven or with an oven using the vial technique: The 768 KF Oven is a conventional tube oven. A weighing spoon is used to transfer the sample to the sample boat, which is automatically transported into the oven at the start of the determination. The 832 KF Thermoprep is a further development of the conventional Karl Fischer drying oven. The samples are no longer transported into the oven by using sample boats, but the vial technique is used. The samples are weighed out directly into vials that are then hermetically sealed and placed in

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the oven. A double hollow needle pierces the vial septum. A stream of dry carrier gas (air or inert gas) is introduced through the inlet needle and passes through the heated sample. The carrier gas containing the released moisture then passes through the outlet needle and via a heated transfer tube directly into the titration cell, in which the Karl Fischer water determination takes place. Depending on the water content of the sample the determination can be carried out volumetrically or, in the trace range, coulometrically. Working with vials means that the oven cannot be contaminated by the sample and, as a result, no carryover or memory effects can occur. In addition, the hermetic vial seals with PTFE-coated septa reliably prevent the water content from being falsified by atmospheric moisture.

6.4 Automation in Karl Fischer titrations The use of a sample changer is well worth while for large numbers of samples. The advantages are not only to be found in the amount of time saved by the laboratory personnel: automated systems control the operating procedures and improve both the reproducibility and accuracy. In automated KF titration there is an additional difficulty; normal KF sample changers do not work in closed systems and penetration by atmospheric moisture cannot be prevented. For a water content >1% this interference can be corrected by the subtraction of a previously determined blank value. In principle, the use of such sample changers can only be recommended for volumetric determinations. The predefined

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methods of the 824 Easy Sample Changer enable the problem-free and simple processing of a series of samples. The comprehensive Karl Fischer equipment allows the determination of water in both solids and liquids. The samples are weighed out into glass vessels and covered with aluminum foil. For the determination the 824 Easy Sample Changer transfers the particular sample to the titration station and moves the titration head to the working position. A remote signal starts the addition of the solvent and then the titration itself. However, Metrohm also offers you the possibility of the automated determination of samples with a water content $<1\%$: The 789 Sample Processor XL is the first sample changer in the world that supplies reliable and reproducible results in combination with the KF Coulometer. As in the oven technique, the samples are weighed out into vials, hermetically sealed with septa and placed on the sample rack. With the aid of a pipetting system adapted to the special conditions of the KF titration (minimization of the penetration of atmospheric humidity) a defined volume of the sample is transported to the coulometer cell where the water content is determined. By using the 774 Oven Sample Processor the automation of the heating method is also possible. The principle corresponds to that of the 832 KF Thermoprep, with the exception that the sample vials are placed on the rack of the Oven Sample Processor. During the analysis the turntable positions the sample vial above the oven. The vial is then

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lowered automatically into the heating block. In addition, the Oven Sample Processor offers the possibility of using temperature gradients. The recorded water release curve can be used to determine the optimum heating temperature for the particular sample. The curve also allows statements to be made about the kinetics of water release as a function of the temperature. Depending on the water content of the sample, the titration can be carried out volumetrically or, in the trace range, coulometrically.

6.5 Accessories for Karl Fischer titrations

A wide range of accessories makes it easier to carry out Karl Fischer titrations in the laboratory: Use the 703 KF Titration Stand for stirring and manual exchange of the used titration solution. Solvents can be added and aspirated off with the built-in membrane pump without having to open the cell. This shortens the conditioning time considerably and minimizes reagent handling by the user. The 803 Ti Stand is similar to the 703 Ti Stand, but its functions and design are adapted to suit Titrados. With KF Coulometers and a Titrando, reagent exchange is even easier to carry out by using a Dosino. All you have to do is define the time and the instrument will carry out the exchange automatically. If special additional steps are necessary, such as titration of a heated solution or homogenization of the sample in the titration cell, then special accessories are available to suit

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your particular application. Some examples:

7. Side reactions

KF water determination is only specific if no side reactions take place with the KF reagents. This means that no water should be released in side reactions, nor should the sample consume or release iodine. Most (unwanted and interfering) side reactions can be suppressed by suitable measures. The most important interfering substances are:

7.1. Carbonates, hydroxides and oxides

During the KF titration methyl sulfuric acid and hydriodic acid are produced. Although these are attached to the base, they still form weak acids that react with the above-mentioned compounds to form salts and release water: Not all oxides take part in an acidic reaction. In addition to alkali oxides and alkaline earth oxides, the following will react: We recommend that the H₂O in these compounds is determined by the oven method.

7.2. Aldehydes and ketones

These two classes of compound combine with the methanol in the KFR to form acetals and ketals respectively with the release of water. With aldehydes the reactivity decreases as the chain length increases. Aromatic aldehydes react more slowly than aliphatic aldehydes. Acetaldehyde reacts the quickest and causes the most problems. But there

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are always exceptions to any rule, formaldehyde and chloral do not form acetals and can be titrated without any problems.

Ketones are less reactive than aldehydes. The reactivity also decreases as the chain length increases. Ketal formation is most rapid with acetone and cyclohexanone. The following are regarded as being stable (no ketal formation): diisopropyl ketone, benzophenone, deoxybenzoin, benzil, benzoin, camphor, alizarin, dibenzalacetone and dichlorotetrafluoroacetone.

Acetal and ketal formation is also influenced by the alcohol present. As the chain length increases the reactivity also decreases. This means that methanol is the quickest to react. There are several ways of preventing acetal or ketal formation: a) By using KF reagents specially labeled for this purpose. b) By using methanol-free KF reagents. Methanol is replaced by a higher alcohol, e.g. 2-methoxyethanol. c) By making use of the different reaction rates and slowing down acetal or ketal formation so much that they do not interfere with the titration. This can be achieved by carrying out the titration at $-10\text{ }^{\circ}\text{C}$. This also freezes the bisulfite addition (see following section). With aldehydes a second reaction – bisulfite addition – can also interfere. In this case H_2O is consumed, so that the water content determined could be too low.

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7.3. Strong acids Strong acids can react with the methanol contained in the KF reagent and release water: Among the organic acids the strongest esterification is shown by formic acid; acetic acid only undergoes partial esterification. Methanol esterifies the quickest. With propanol and/ or 2-methoxyethanol slower reaction rates are achieved. However, strong acids also lower the pH of the reaction solution so much that the KF reaction either no longer takes place or proceeds unsatisfactorily. The reaction mixture can be buffered by the addition of a base (e.g. imidazole) and the pH brought back to the optimal range again. If this measure is taken then no esterification occurs even with strong mineral acids

7.4. Silanols/Siloxanes The terminal silanol groups esterify with the methanol contained in the KF reagents and release water. This means that methanol-free reagents must be used to determine the water content of these substances

7.5. Boron compounds Boric acid reacts with methanol with esterification and the release of water: Boron trioxide (B_2O_3), metaboric acid (HBO_2) and borax / sodium tetraborate ($Na_2B_4O_7$) all react in the same way. The esterification rate decreases as the chain length of the alcohol increases (methanol esterifies the quickest).

7.6. Metal peroxides These form water with KF solutions and cannot be titrated directly. In some cases the oven method can be used (provided

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that the substances release their H₂O at <250 °C).

7.7. Reducing agents Reducing agents can react with the iodine in the KF reagents and simulate a water content that is too high. 0.5 mol H₂O / mol substance in the sample is simulated by

Petrochemical products may also contain interfering substances. The additive tetraethyl lead simulates a higher water content: The water content of salts can usually be determined by the oven method. We recommend the following procedure for determining the water content of ascorbic acid and hydrazines: The sample is dissolved in methanol and titrated against methanolic iodine solution. This oxidizes it and the KF water determination is no longer interfered with. In the same solution a titration can then be carried out with a combined KF solution. The H₂O content of the methanol used must be determined separately and subtracted from the result.

7.8. (Strong) Oxidizing agents These release elemental iodine from the iodides in the KF reaction and simulate a water content that is too low.

7.8 a Dialkyl peroxides, diacyl peroxides and peresters form iodine. However, the reaction takes place slowly and can be «frozen» at low temperatures. The titration temperature to be used is product-specific and can be as low as -60 °C. Hydroperoxides – e.g. H₂O₂ – do form iodine, but also an equivalent

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amount of H₂O. The KF titration is free from interference

7.8 b Elemental chlorine and bromine also release iodine. In this case excess SO₂ is passed through the sample. This reduces these substances to chloride and bromide respectively, which no longer interfere.

7.8 c The water content of dichromates, chromates and permanganates is determined by the oven method.

8. Sample preparation techniques

Each sample preparation is preceded by sampling. Particularly with inhomogeneous samples – and there are a lot of them, even if it is not always obvious – this initial step plays a crucial role for the correctness and accuracy of the results. Correct sampling is therefore an essential requirement for the subsequent determination! Many standardized methods provide detailed instructions concerning this subject. This may also apply for in-house instructions. If no instructions are available then the following points must be observed:

- The water content of the sample taken must be the same as the average water content of the original sample – i.e. the sample must be representative.
- Atmospheric moisture is an ever-present source of contamination, and not only in tropical countries.
- Hygroscopic samples have a large amount of water adhering to their surfaces.
- Samples that easily release water or become weathered have only a small

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amount of water on their surfaces.

- Liquids can release water if the water is not present in dissolved form but as a dispersion (precipitation or flotation, adherence to vessel walls – e.g. petrochemical products, margarine). The same applies for samples in which the solubility limit for water is undercut on cooling down.
- The more heterogeneous the water distribution in the sample, the larger the sample to be taken.
- Sampling must take place as quickly as possible and the water content of the sample must not change during transport and storage.

General Many samples can be weighed out and/or injected directly into the titration vessel or titration cell, provided that they are soluble in the solvent used and do not undergo side reactions with the KFR. In contrast, other samples require one or more preparation steps. The most important sample preparation techniques are described below. Further information can be found under the individual classes of substances. It is very important that the sample preparation technique used does not introduce any additional water into the sample and that no water losses occur by heating the sample too strongly.

A) Dilution For reasons of accuracy (sample size) we recommend that samples with a high water content (>40%) are diluted with an inert dry solvent before the determination. Methanol can normally be used for this purpose. If the sample contains emulsified fats and/or oils then a dry mixture of 60% methanol and 40%

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decanol (or, if permitted, chloroform) should be used. The procedure is as follows: Approx. 1 g sample is weighed out exactly into a dry Erlenmeyer flask fitted with a septum stopper. This is treated with approx. 20...25 g dry solvent (mixture), the flask is then sealed and the contents mixed. Part of the mixture is drawn off in a dry syringe and the syringe and its contents are tared. Some of the mixture is injected into the titration vessel and titrated. The amount of sample injected is determined by back-weighing the syringe. If a dry solvent is used then it is not absolutely necessary to take its blank value into account. Here is an example:

The following samples can be treated in this way: – Drilling and cutting oil emulsions – Cosmetic lotions and emulsions – Spirits – Milk and cream, ready-to-use salad sauces, liquid spices. B) Solvent extraction

This can take place in two different ways – external or during the titration in the titration vessel. It is usually carried out either warm or hot, less frequently at room temperature. If an external extraction is carried out then please remember that all the water can never be extracted, an average water content between the sample and the solvent is obtained.

Example: 10 g sample with a water content of 10% is extracted with 50 mL methanol (density approx. 08 g/cm³) with a water content of 100 ppm. The H₂O content of the mixture is 1004 mg / 50 g, corresponding to 20.08 mg H₂O / g or 2.008% H₂O. This also means that the larger the ratio between the sample and dry

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solvent, the less water will be extracted from the sample. The blank value of the solvent must, of course, be taken into consideration when calculating the water content. Better values are obtained if the extraction is carried out directly in the titration vessel and the water is titrated immediately.

A titration vessel equipped with a heating jacket and a reflux condenser is used and the titration is carried out, for example, at 50 °C. Work can be carried out with or without an extraction time (set at the titrator) as required. (If an extraction time is used then please note that a sufficiently long stop delay must be observed for the EP.) This method is primarily used for cereal flour. C) Homogenization In liquids water can float, adhere to the inner walls of the sample bottle or be deposited on the bottom. This occurs particularly when the sample is taken at a high temperature and allowed to cool down to room temperature before the analysis. The water can again be distributed homogeneously by immersing the bottle in an ultrasonic bath. The separated water can also be dissolved again by the addition of an exactly weighed amount of a solubility promoter (e.g. isopropanol), whose water content has been accurately determined (blank value).

This method is suitable for pastes and viscous liquids. Solids are usually inhomogeneous and must be thoroughly comminuted and homogenized. The procedure depends on the constitution of the sample:

Coarse particles, hard and tough Ball mill or analytical mill (closed) with water cooling. This method is suitable

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for coffee, cereals, semolina, pasta, beans, lentils, pills, spices* and tea*. Coarse, fatty Grating and high-frequency comminution, e.g. for nuts, seeds and almonds.

Coarse, soft, inhomogeneous Comminute with a meat grinder or knife, then high-frequency comminution. Suitable for dried fruits, meat and sausages, fruit and vegetables.

Pasty, inhomogeneous Homogenization with a blender, then possibly high-frequency comminution. For jam, mustard, precooked spinach.

Suspensions High-frequency comminution, e.g. for fruit and vegetable juice, fruit pulp, concentrates and tomato puree. High-frequency comminution is the only method that can guarantee that even plant samples are thoroughly homogenized. In this method the cells are disrupted and the water present inside the cells is extracted. Comminution is usually carried out externally in a dry Erlenmeyer or round-bottom flask with the addition of a suitable solvent (mostly methanol). Comminution in the titration vessel is also possible; the mixing head is commonly used as a stirrer (remove the magnetic stirrer bar from the titration vessel).

D) Heating method The sample is either heated in a dry stream of gas (air or nitrogen) or a dry stream of gas is passed through the heated sample; in each case the released water is transferred to the titration vessel and immediately titrated. The heating

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temperature should be set as high as possible – just high enough so that the sample does not decompose and no interfering secondary constituents distill over. The heating temperature and duration must be optimized in preliminary tests. The appropriate blank values must be subtracted. This method is used for samples that cannot be extracted (or cannot be completely extracted) with a suitable solvent or that undergo side reactions with the KFR. Examples are: – Plastics, minerals, carbonates and metal oxides, carbon powder – Petrochemical products containing additives (motor oils, lubricating oils and greases) – Starting materials for pharmaceutical products, certain freeze-dried samples.

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